

Monthly radioactivity measurement "June, 2012" was conducted.

Measurement Result:

No radioactive substances (Iodine/Caesium) was detected with any of Products/Cannulae.

Table 1. Radioactivity measurement - June/2012

Sample Category	Nuclide				Radiation Dose	
	Iodine	Caesium				
Unit	I-131	Cs-134	Cs-137	Cs-136	μSv/h	cpm
Products (Needle)	ND ^{*1}	ND ^{*1}	ND ^{*1}	ND ^{*1}	—	—
Cleaning Water	Measurement once in two months (Next July)					
Cannulae	Measurement once in two months (Next August)					
Components	—	—	—	—	ND ^{*1}	—

*1) Not Detected (Equivalent count rate as Background)